



## **SonicWeld** Rx<sup>®</sup>

The Perfect Choice



Oral and maxillo-facial surgery is our passion! Its further development, together with our customers, is our ambition. Every day we work on developing innovative products and services which meet the highest demands on quality, and which contribute to the wellbeing of the patient.

## Table of Contents

	Pages
Product Features	4-5
Feature, Function, Benefit	6-15
Indications	16-17
References	18-19



## SonicWeld Rx®

It was back in 2001 when KLS Martin launched the osteosynthesis system Resorb x®. Thus offering the first completely resorbable implants made of pure PDLLA. But this was just the beginning.

In 2005, KLS Martin proceeded to revolutionize the field of resorbable osteosynthesis by introducing SonicWeld Rx®, the unique ultrasound technology for insertion of SonicPins.

In 2013, a new chapter in the company's history of resorbables was opened by the introduction of Resorb xG, a PLLA-PGA polymer with improved mechanical features.

Now, KLS Martin is setting up another milestone: The second generation of SonicWeld Rx®. The novel device is an optical highlight in every OR, offering improved and additional features for a user-friendly application. Just see for yourself.

**SonicWeld Rx®. The perfect choice.**

**2015:** SonicWeld Rx®  
second generation


**2013:** Resorb xG

**2005:** SonicWeld Rx®

**2001:** Resorb x®



## Feature, Function and Benefit



SonicWeld Rx® is a revolutionary technique for use in craniomaxillofacial osteosynthesis. It combines highly advanced ultrasound technology with resorbable implants to provide extremely stable fixation and completely eliminate the need for a second operation.

The procedure is simple: resorbable meshes are heated up, shaped to fit the application site and then fixed in place with SonicPins inserted into predrilled holes. This is done with a sonotrode that liquefies the pins, thus causing them to bond with the meshes and penetrate into the bone cavities to anchor themselves securely.

The method is clinically certified and validated and very patient-friendly as well. The implants degrade through natural hydrolysis in a controlled process. SonicWeld Rx® is primarily stable, convenient, fast, easy and safe. Designed for cranial fixation, ideal for pediatric trauma, and indicated also for cancellous bone structures.

	Feature and Function	Benefit
	<ul style="list-style-type: none"><li>■ The ultrasonic energy sets the SonicPin into mechanical vibration</li><li>■ The liquid SonicPin penetrates into the bone cavities</li><li>■ The principle works both in cortical and cancellous spongy bone</li><li>■ Low power effort during SonicPin insertion</li><li>■ Implantation of the SonicPin in angle position is possible</li><li>■ Maximum temperature increase of the bone at about 1 mm from the implant: 11 °C</li><li>■ Only 30 - 40 seconds after SonicPin insertion, temperature increase is below 5 °C</li><li>■ No risk of pin/screw breakage</li></ul>	<ul style="list-style-type: none"><li>■ The material liquifies at the interface between the pre-drilled bone and the SonicPin via friction</li><li>■ The material reaches bone cavities beyond the reach of common screws</li><li>■ Excellent three-dimensional stability both in cortical and spongy bone</li><li>■ Particularly effective in poorer bone quality</li><li>■ Repositioning of small bone fragments</li><li>■ Especially suitable in cramped corners without dislocation</li><li>■ Maximum bone temperature is below denaturing temperature of 56 °C</li><li>■ No bone necrosis</li><li>■ Fast cooling down of the material and surrounding bone</li><li>■ Secure anchorage of the SonicPin in the bone only three seconds after activation</li><li>■ No emergency system is necessary</li></ul>
 <p>in vivo Tensile Strength</p> <p>in vivo Side Force Bending</p> <p>■ Resorbable screws ■ SonicPin Rx</p>	<ul style="list-style-type: none"><li>■ Locking effect between the SonicPin and the pre-drilled hole</li><li>■ Locking effect between the SonicPin head and the plate</li><li>■ Locking mechanism can be reversed by drilling through the inserted SonicPin</li></ul>	<ul style="list-style-type: none"><li>■ Due to the double locking mechanism extremely stable fixation of the SonicPin in the pre-drilled hole</li><li>■ With SonicPins twice the strength compared to resorbable screws can be achieved</li><li>■ Simple implant removal</li><li>■ Simple correction of the implant position</li></ul>
 <p>Time required to place 20 screws/pins</p> <p>■ Resorbable screws ■ SonicPin Rx</p>	<ul style="list-style-type: none"><li>■ No need for pre-tapping</li></ul>	<ul style="list-style-type: none"><li>■ Exceptionally fast implantation of the SonicPin</li><li>■ Reduction in surgical time</li></ul>

## Feature, Function and Benefit



Two resorbable polymers for osteosynthesis, PDLLA and PLLA-PGA, have been well-established in craniomaxillofacial surgery.

**Resorb x®** polymer is a 100% Poly-D,L-Lactic Acid (PDLLA).

**Resorb xG** polymer consists of 85% Poly-L-Lactic Acid (PLLA) and 15% Poly Glycolic Acid (PGA).

Both resorbables maintain the majority of their strength for 8-10 weeks, allowing complete fracture healing and bone regeneration.

The core of the degradation process:

The complex polymer chains absorb the water contents (H<sub>2</sub>O molecules) of surrounding body fluids through a process called “hydrolysis”.


The stored water initiates the degradation process by continuously breaking down the long polymer chains into ever shorter structures or simpler molecules. Metabolic pathways subsequently transform the molecules into carbon dioxide and water; both of these compounds are discharged naturally.

## SonicWeld Rx®

### Resorb x®


#### Feature and Function

#### Benefit




- Polymer consists of 100% Poly-D,L-Lactic Acid (PDLLA)

- Totally amorphous polymer
- Residue free degradation
- Numerous animal and clinical studies prove excellent biocompatibility and a safe degradation process.
- Resorption time observed in ultrasound follow-up: 12 - 30 months




### Resorb xG




- Polymer consists of 85% Poly-L-Lactic Acid (PLLA) and 15% Poly Glycolic Acid (PGA)

- Higher initial strength
- Faster decrease of both strength and mass
- Resorption time: approximately 12 - 14 months



## Feature, Function and Benefit




SonicPins are characterized by their unique geometry. The geometry guarantees maximum polymer outflow in the surrounding bone cavities during SonicPin insertion. Thus reducing the power effort for SonicPin insertion to a minimum. Sonic Pins are available in two diameters:

- **green clip:** Ø 1.6 mm
- **red clip:** Ø 2.1 mm

Resorbable implants are available in various designs and thicknesses to give the surgeon options to match every indication. The holes of the plates and meshes are perfectly adapted to the geometry of the SonicPins. Thus the head of the SonicPin is optimally countersunk in the implant.

## SonicWeld Rx®

### SonicPins




#### Feature and Function

- Color-coded clip magazines
  - green: SonicPins Ø 1.6 mm
  - red: SonicPins Ø 2.1 mm
- Self-retaining pin head
- Optimized pin geometry
- Both SonicPin sizes fit all implants of Resorb x® and Resorb xG product range
- Sterile delivery

#### Benefit

- Easy identification of the appropriate SonicPin diameter
- Convenient pin removal from clip magazine
- Maximum polymer outflow in the surrounding bone structure
- Easy pin insertion
- Complete cross compatibility
- Always ready to use



### SonicPin types



- Standard SonicPin
- Micro SonicPins Rx without pin head
- Endobrow SonicPins Rx with specially designed pin tip for sutures

- Perfect solution for a wide range of applications
- Ideal for narrow spaces, e. g. preprosthetic augmentation
- Ideal for endobrow lifting

### Plates, Meshes, Foils and Membranes



- Huge variety of different geometries, sizes and thicknesses
- Round edge geometry
- Can easily be contoured in the Xcelsior water bath and cut with scissors intraoperatively
- Flexible meshes
- Membranes and foils with minimal thickness (0.1, 0.2 or 0.3 mm)
- All Resorb x® and Resorb xG implants fit both SonicPin diameters (1.6 and 2.1 mm)
- Sterile delivery

- Right implant for every indication
- Minimal palpability and susceptibility
- Easy adaption to patient-specific anatomy
- Very easy to adapt to patient specific anatomy
- Ideal for preprosthetic augmentation
- Complete cross compatibility
- Always ready to use

## Feature, Function and Benefit




The ultrasonic unit of the SonicWeld Rx® system converts electric energy into mechanical vibrations (ultrasound).

When using a standard sonotrode, the ultrasonic energy causes a phase change of the resorbable material at the interfaces between the bone and the SonicPins via friction. Thus the SonicPin glides into the predrilled hole. When using a smoothing sonotrode, the ultrasonic energy allows to smooth the resorbable implants (e. g. a membrane).

## SonicWeld Rx®

### Ultrasonic unit



#### Feature and Function

- Simple and elegant design
- Round edge geometry
- Two handles to carry the device
- Two connecting sockets for handpieces
- One pre-defined power level
- Opportunity to choose the individual system language

#### Benefit

- Clear optical distinction to first generation device
- Easy to clean
- Secure fit of the device during transportation
- Possibility to work alternately with two sonotrodes (e.g. a standard and a smoothing sonotrode or two standard sonotrodes)
- Optimal system setting
- User-friendly application
- No comprehensive problems


### Handpiece



- Ergonomically designed handpiece
- Finger activation
- Light and acoustic support during activation
- Autoclavable

- Well balanced and comfortable fit
- Exclusive concentration on the hand during SonicPin insertion or smoothing
- 1:1 feedback during activation period
- Guaranteed biocompatibility for 250 sterilization cycles

### Sonotrodes



- Standard sonotrodes
  - straight
  - angled
- Smoothing sonotrodes
  - straight
  - angled

- Ideal for SonicPin insertion in straight position
- Combined sonotrode  
Ideal for SonicPin insertion in angled position (e. g. orbita or side tooth area)
- Smoothing of implants in straight position
- Smoothing of implants in straight or angled position (e. g. orbita or side tooth area)

## Feature, Function and Benefit




The Xcelsior water bath is intended for heating up resorbable implants for the purpose of adapting them to the patient's anatomical conditions (e. g. bone surface). Various templates are available that help to adapt the implants to the shape of the bone.

The BOS drill is a fully-fledged and universally applicable drill system. The battery tools do not require a charger or base unit and are always ready – wherever and whenever you need them.

## SonicWeld Rx®

### Xcelsior water bath




#### Feature and Function

- Tool for heating up Resorb x® and Resorb xG implants in the hot water (70 - 90 °C / 158 - 194 °F) to adapt it to the patient-specific bone contour
- Sterilizable material ① ②

#### Benefit

- Perfect temperature range to adapt Resorb x® and Resorb xG implants
- To be used in the sterile area of the OR


### Templates



- Various templates available
- Adaption of the implant to the patient's anatomical condition in the Xcelsior water bath

- Template reflects the implant 1-to-1
- Safe selection of the sterile-packed implant
- Perfect fit of implant


### BOS Drill



- 600 rev/min, high-speed forward
- Ergonomic design
- Lightweight handle weighing only 200 g
- Can be operated with a finger
- Sterile battery pack simply needs to be clicked-on

- Ideal for predrilling
- Safe fit in the user's hand
- Especially indispensable when dealing with a large number of implants
- Comfortable to use
- Always charged and ready for use

## Step by Step to innovative Osteosynthesis





### Indications

The KLS Martin Resorb x® and Resorb xG implants are intended for surgical procedures in which an internal fixation by resorbable implants is required for aligning, reconstructing and stabilizing bone tissue.



Craniofacial corrective  
osteotomies  
(e. g. craniosynostosis)



Osteosynthesis in non-load-bearing areas of the craniomaxillofacial skeleton

Preprosthesis augmentation

Endobrow fixation

## Literature

Abdel-Galil, K. & Loukota, R.,  
Fixation of comminuted diacapitular fractures  
of the mandibular condyle with ultrasoand-  
activated resorbable pins.  
Br J Oral Maxillofac Surg, 46(6), 2008,  
S. 482-484

Aldana, P. R., et al.,  
Ultrasoand-aided fixation of a biodegradable  
cranial fixation system: uses in pediatric  
neurosurgery.  
J Neurosurg Pediatr, 3(5), 2009, S. 420-424

Aldana, P. R., et al.,  
Ultrasoand-aided fixation of biodegradable  
implants in pediatric craniofacial surgery.  
Pediatr Neurosurg, 47(5), 2011, S. 349-353

Arnaud, E. & Renier, D.,  
Pediatric craniofacial osteosynthesis and  
distraction using an ultrasonic-assisted pinned  
resorbable system: a prospective report with  
a minimum 30 months' follow-up.  
J Craniofac Surg, 20(6), 2009, S. 2081-2086

Basa, S.,  
Does ultrasonic resorbable pin fixation offer  
predictable results for augmentation eminoplasty  
in recurrent dislocations?  
J Oral Maxillofac Surg, 72(8), 2014,  
S. 1468-1474

Buijs, G. J., et al.,  
Mechanical strength and stiffness of the bio-  
degradable SonicWeld Rx® osteofixation system.  
J Oral Maxillofac Surg, 67(4), 2009, S. 782-787

Burger B. W.,  
Use of ultrasoand-activated resorbable poly-D,  
L-lactide pins (SonicPins) and foil panels  
(Resorb x®) for horizontal bone augmentation  
of the maxillary and mandibular alveolar ridges.  
J Oral Maxillofac Surg, 68(7), 2010,  
S. 1656-1661

Chen, Y. B. & Zhang, H. Z.,  
Ultrasoand-aided biodegradable osteosynthesis  
system: application in fixation of oral and  
maxillofacial fractures.  
Zhonghua Yi Xue Za Zhi, 93(18), 2013,  
S. 1418-1421

Cho P. W. J., et al.,  
Biomechanical study of SonicWeld Rx® pin  
in cortical bone graft layering technique.  
J Oral Maxillofac Surg, 69(5), 2011,  
S. 1519-1524

Cristofaro, M. G., et al.,  
A new system of resorbable rigid three-  
dimensional fixation using ultrasoands  
(SonicWeld Rx®+ Sonic Pins Rx) adopted  
in craniofacial traumatology: the author's  
experience.  
It J Maxillofac Surg, 20, 2009, S. 4-52

Eckelt U., et al.,  
Ultrasoand aided pin fixation of biodegradable  
osteosynthetic materials in cranioplasty for  
infants with craniosynostosis.  
J Craniomaxillofac Surg, 35(4-5), 2007,  
S. 218-221

Heidemann W., et al.,  
Degradation of poly(D,L)lactide implants with  
or without addition of calciumphosphates in vivo.  
Biomaterials, 22(17), 2001, S. 2371-2381

Heidemann W. & Gerlach K. L.,  
Sonographic examinations on the degradation  
of bioresorbable osteosynthesis materials.  
Biomed Tech, 46(9), 2001, S. 236-240

Heidemann W, et al.,  
In vivo investigation of the degradation of  
poly(D,L)lactide and poly(L-lactide-co-glycolide)  
osteosynthesis material.  
Mand-, Kiefer- GesichtsChir, 7, 2003,  
S. 283-288

Iglhaut, G.,  
The Minimally Invasive Shell Technique  
for Bone Augmentation.  
Oralchirurgie Journal, 9(3), 2009

Iglhaut, G., et al.,  
Shell technique using a rigid resorbable barrier  
system for localized alveolar ridge augmentation.  
Clin Oral Implants Res, 25(2), 2014,  
S. 149 - 154

Lee, J. H. & Park, J. H.,  
The clinical usefulness of ultrasound-aided fixation  
using an absorbable plate system in patients with  
zygomatico-maxillary fracture.  
Arch Plast Surg, 40(4), 2013, S. 330-334

Mai, R.,  
Bone welding – a histological evaluation in the jaw.  
Ann Anat, 189(4), 2007, S. 350-355

Meara, D. J., et al.,  
Fixation of Le Fort I osteotomies with poly-DL-lactic  
acid mesh and ultrasonic welding –  
a new technique.  
J Oral Maxillofac Surg, 70(5), 2012, S. 1139-1144

Müller-Richter, U. D., et al.,  
Treatment of intracapsular condylar fractures  
with resorbable pins.  
J Oral Maxillofac Surg, 69(12), 2011,  
S. 3019-3025

Pietrzak W.S.,  
Bioabsorbable polymer applications in  
musculoskeletal fixation and healing.  
In: Pietrzak W. S. (ed.), Orthopedic biology  
and medicine: Musculoskeletal tissue  
regeneration, biological materials and  
methods,  
Totawa: Humana Press, 2008, S. 509-529

Pilling E., et al.,  
An Experimental study of the biomechanical  
stability of ultrasound-activated pinned  
(SonicWeld Rx® + Resorb x®) and screwed  
fixed (Resorb x®) resorbable materials for  
osteosynthesis in the treatment of simulated  
craniosynostosis in sheep.  
Br J Oral Maxillofac Surg, 45(6), 2007,  
S. 451-456

Pilling, E., et al.,  
An experimental in vivo analysis of the  
resorption to ultrasound-activated pins  
(SonicWeld Rx®) and standard biodegradable  
screws (Resorb x®) in sheep.  
Br J Oral Maxillofac Surg, 45(6), 2007,  
S. 447-450

Reichwein, A.,  
Clinical experiences with resorbable  
ultrasonic-guided, angle-stable osteosynthesis  
in the panfacial region.  
J Oral Maxillofac Surg, 67(6), 2009,  
S. 1211-1217

Schneider, M., et al.,  
Stability of fixation of diacapitular  
fractures of the mandibular condylar  
process by ultrasound-aided resorbable pins  
(SonicWeld Rx® System) in pigs.  
Br J Oral Maxillofac Surg, 49(4), S. 297-301

Schneider, M., et al.,  
Ultrasound-aided resorbable osteosynthesis  
of fractures of the mandibular condylar base:  
an experimental study in sheep.  
Br J Oral Maxillofac Surg, 50(6), 2012,  
S. 528-532

Stelnicki, E. J., et al.,  
Use of absorbable poly (D,L) lactic acid plates  
in cranial-vault remodeling: presentation of the  
first case and lessons learned about its use.  
Cleft Palate Craniofac J, 42(4), 2005,  
S. 333-339

Völker, W., et al.,  
The use of resorbable osteosynthesis materials.  
Laryngorhinootologie. 90(1), 2011, S. 23-25

Wood, R. J., et al.,  
New resorbable plate and screw system  
in pediatric craniofacial surgery.  
J Craniofac Surg, 23(3), 2012, S. 845-849

## KLS Martin Group

### KLS Martin Australia Pty Ltd.

Sydney · Australia  
Tel. +61 2 9439 5316  
australia@klsmartin.com

### KLS Martin do Brasil Ltda.

São Paulo · Brazil  
Tel. +55 11 3554 2299  
brazil@klsmartin.com

### KLS Martin Medical (Shanghai) International Trading Co. Ltd.

Shanghai · China  
Tel. +86 21 5820 6251  
china@klsmartin.com

### KLS Martin India Pvt Ltd.

Chennai · India  
Tel. +91 44 66 442 300  
india@klsmartin.com

### Martin Italia S.r.l.

Milan · Italy  
Tel. +39 039 605 67 31  
italia@klsmartin.com

### Nippon Martin K.K.

Tokyo · Japan  
Tel. +81 3 3814 1431  
nippon@klsmartin.com

### KLS Martin SE Asia Sdn. Bhd.

Penang · Malaysia  
Tel. +604 505 7838  
malaysia@klsmartin.com

### KLS Martin de México S.A. de C.V.

Mexico City · Mexico  
mexico@klsmartin.com

### Martin Nederland/Marned B.V.

Huizen · The Netherlands  
Tel. +31 35 523 45 38  
nederland@klsmartin.com

### Gebrüder Martin GmbH & Co. KG

Moscow · Russia  
Tel. +7 499 792-76-19  
russia@klsmartin.com

### KLS Martin Taiwan Ltd.

Taipei 106 · Taiwan  
Tel. +886 2 2325 3169  
taiwan@klsmartin.com

### Gebrüder Martin GmbH & Co. KG

Dubai · United Arab Emirates  
Tel. +971 4 454 16 55  
middleeast@klsmartin.com

### KLS Martin UK Ltd.

London · United Kingdom  
Tel. +44 1189 000 570  
uk@klsmartin.com

### KLS Martin LP

Jacksonville · Florida, USA  
Tel. +1 904 641 77 46  
usa@klsmartin.com

### Gebrüder Martin GmbH & Co. KG

#### A company of the KLS Martin Group

KLS Martin Platz 1 · 78532 Tuttlingen · Germany  
P.O. Box 60 · 78501 Tuttlingen · Germany  
Tel. +49 7461 706-0 · Fax +49 7461 706-193  
info@klsmartin.com · www.klsmartin.com




06.19 · 90-421-02-08 · Printed in Germany · Copyright by Gebrüder Martin GmbH & Co. KG · Alle Rechte vorbehalten · Technische Änderungen vorbehalten  
We reserve the right to make alterations · Cambios técnicos reservados · Sous réserve de modifications techniques · Ci riserviamo il diritto di modifiche tecniche



The KLS Martin SonicWeld Rx<sup>®</sup> solution is based on the BoneWelding<sup>®</sup> technology<sup>1)</sup> protected by the industrial property rights of WoodWelding AG, Switzerland, and has been licensed by this company.

<sup>1)</sup> "BoneWelding" is a registered Swiss trademark.



## **SAM**

Standard Approach Mesh  
Pre-contoured titanium mesh implants



Oral and maxillo-facial surgery is our passion! Its further development, together with our customers, is our ambition. Every day we work on developing innovative products and services which meet the highest demands on quality, and which contribute to the wellbeing of the patient.


## **SAM** – Standard Approach Mesh Pre-contoured titanium mesh implants

### Benefits


- Pre-contoured titanium mesh implants and templates for 11 standard neurosurgical approaches
- Mesh sizes and shapes were created using a compilation of adult CT scans to develop the best possible “average” fit
- 0.6 mm thick titanium: rigid, yet able to modify slightly to match each patient’s specific anatomy
- Used in conjunction with KLS Martin Drill-Free 1.5 mm Neuro maxDrive® screws
- Each implant comes with detachable identification tag to aid in orientation and to maximize inventory control (see below)
- Also available in sterile version









Blue color designates frontal direction of mesh




25-007-05-09  
contoured mesh  
with I.D. tag attached
















25-007-01-09    
 25-007-01-71    
 Mesh for bifrontal approach  
 = 0.6 mm


## **SAM** – Standard Approach Mesh Pre-contoured titanium mesh implants









25-007-02-09    
 25-007-02-71    
 Mesh for pterional approach, left  
 = 0.6 mm









25-007-04-09    
 25-007-04-71    
 Mesh for parietal approach, left  
 = 0.6 mm









25-007-06-09    
 25-007-06-71    
 Mesh for anterior parasagittal approach  
 = 0.6 mm









25-007-03-09    
 25-007-03-71    
 Mesh for pterional approach, right  
 = 0.6 mm



25-007-05-09    
 25-007-05-71    
 Mesh for parietal approach, right  
 = 0.6 mm



25-007-07-09    
 25-007-07-71    
 Mesh for posterior parasagittal approach  
 = 0.6 mm



25-007-00-09

Mesh set, pre-contoured, complete

= 0.6 mm


25-007-70-09

Template set

Explanation of icons

- Pure titanium
- Packaging unit
- Plate profile

**STERILE | R** Sterile packed implants




25-007-08-09

25-007-08-71

Mesh for suboccipital approach, left

= 0.6 mm




25-007-10-09

25-007-10-71

Mesh for frontotemporoparietal approach, left

= 0.6 mm




25-007-13-09

25-007-13-71

Mesh for frontotemporoparietal approach, left, XL

= 0.6 mm




25-007-09-09

25-007-09-71

Mesh for suboccipital approach, right

= 0.6 mm




25-007-11-09

25-007-11-71

Mesh for frontotemporoparietal approach, right

= 0.6 mm



25-007-14-09

25-007-14-71

Mesh for frontotemporoparietal approach, right, XL

= 0.6 mm

## **SAM** – Standard Approach Mesh Pre-contoured titanium mesh implants



55-969-32-04 Mesh tray 30 x 30 cm with lid and silicone mat

### Standard Approach Mesh

25-007-00-09	Mesh set, pre-contoured, 0.6 mm complete, consisting of:	Qty/unit
25-007-01-09	Mesh for bifrontal approach	1
25-007-02-09	Mesh for pterional approach, left	1
25-007-03-09	Mesh for pterional approach, right	1
25-007-04-09	Mesh for parietal approach, left	1
25-007-05-09	Mesh for parietal approach, right	1
25-007-06-09	Mesh for anterior parasagittal approach	1
25-007-07-09	Mesh for posterior parasagittal approach	1
25-007-08-09	Mesh for suboccipital approach, left	1
25-007-09-09	Mesh for suboccipital approach, right	1
25-007-10-09	Mesh for frontotemporoparietal approach, left	1
25-007-11-09	Mesh for frontotemporoparietal approach, right	1

All meshes are available in sterile version



55-969-31-04 Mesh tray 30 x 60 cm with lid and silicone mat

## Templates

25-007-70-09	Mesh Template Set complete, consisting of:	Qty/unit
25-007-51-09	Mesh Template for bifrontal approach	1
25-007-52-09	Mesh Template for pterional approach, left	1
25-007-53-09	Mesh Template for pterional approach, right	1
25-007-54-09	Mesh Template for parietal approach, left	1
25-007-55-09	Mesh Template for parietal approach, right	1
25-007-56-09	Mesh Template for anterior parasagittal approach	1
25-007-57-09	Mesh Template for posterior parasagittal approach	1
25-007-58-09	Mesh Template for suboccipital approach, left	1
25-007-59-09	Mesh Template for suboccipital approach, right	1
25-007-60-09	Mesh Template for frontotemporoparietal approach, left	1
25-007-61-09	Mesh Template for frontotemporoparietal approach, right	1

## KLS Martin Group

### **KLS Martin Australia Pty Ltd.**

Sydney · Australia  
Tel. +61 2 9439 5316  
australia@klsmartin.com

### **KLS Martin do Brasil Ltda.**

São Paulo · Brazil  
Tel. +55 11 3554 2299  
brazil@klsmartin.com

### **KLS Martin Medical (Shanghai) International Trading Co., Ltd.**

Shanghai · China  
Tel. +86 21 5820 6251  
china@klsmartin.com

### **KLS Martin India Pvt Ltd.**

Chennai · India  
Tel. +91 44 66 442 300  
india@klsmartin.com

### **Martin Italia S.r.l.**

Milan · Italy  
Tel. +39 039 605 67 31  
italia@klsmartin.com

### **Nippon Martin K.K.**

Tokyo · Japan  
Tel. +81 3 3814 1431  
nippon@klsmartin.com

### **KLS Martin SE Asia Sdn. Bhd.**

Penang · Malaysia  
Tel. +604 505 7838  
malaysia@klsmartin.com

### **Martin Nederland/Marned B.V.**

Huizen · The Netherlands  
Tel. +31 35 523 45 38  
nederland@klsmartin.com

### **Gebrüder Martin GmbH & Co. KG**

Moscow · Russia  
Tel. +7 499 792-76-19  
russia@klsmartin.com

### **Gebrüder Martin GmbH & Co. KG**

Dubai · United Arab Emirates  
Tel. +971 4 454 16 55  
middleeast@klsmartin.com

### **KLS Martin UK Ltd.**

London · United Kingdom  
Tel. +44 1189 000 570  
uk@klsmartin.com

### **KLS Martin LP**

Jacksonville · Florida, USA  
Tel. +1 904 641 77 46  
usa@klsmartin.com

### **Gebrüder Martin GmbH & Co. KG**

#### **A company of the KLS Martin Group**

KLS Martin Platz 1 · 78532 Tuttlingen · Germany  
P.O. Box 60 · 78501 Tuttlingen · Germany  
Tel. +49 7461 706-0 · Fax +49 7461 706-193  
info@klsmartin.com · www.klsmartin.com

